

# BIODECONTA

## 去污生物指示剂

### 生物指示剂

型号: BDBI-V6      规格:  $10^6$ cfu/pcs  
包装: 50 片 / 包      保质期: 9 个月  
菌种: 嗜热脂肪地芽孢杆菌  
(*Geobacillus stearothermophilus* ATCC 12980)  
D 值: 1.8min (25°C, 500±30 ppm, 70%±1 RS)  
特点: 独立包装, 使用方便; 载体为医用级不锈钢片;  
Tyveke 包装满足过氧化氢穿透要求, 同时避免外部污染。



### 配套培养基

型号: BDBIM01      规格: 1.6ml/ 支  
包装: 50 支 / 盒      颜色: 紫色  
形态: 液体      保质期: 9 个月  
特点: 改良 TSB 培养基, 含 pH 指示剂, 颜色变化明显,  
培养 48 小时可出结果。



本产品专用于汽化过氧化氢 (VHP) 生物去污染过程的监测, 适用于以下场景:

- △ 隔离器、无菌物料传递设备等密闭空间的生物去污效果验证
- △ 洁净厂房、实验室等环境的消毒效果评估



BIODECONTA

# BIODECONTA

## Decontamination biological indicator

### Biological indicators

**Model:** BDBI-V6

**Specifications:**  $10^6$ cfu/pcs

**Packaging:** 50 pieces/pack

**Shelf life:** 9 months

**Bacteria:** *Bacillus dehydrogen*

(*Geobacillus stearothermophilus* ATCC 12980)

**D value:** 1.8min (25°C, 500±30 ppm, 70%±1 RS)

**Features:** Independent packaging, easy to use; the carrier is medical grade stainless steel sheet; avoid external contamination.



### Supporting culture medium

**Model:** BDBIM01

**Specification:** 1.6ml/ pcs

**Packaging:** 50 pieces/pack

**Shelf life:** 9 months

**Form:** Liquid

**Features:** Improved TSB culture medium, containing pH indicator, and the color changes are obvious.

The results will be produced after 48 hours of incubation.



**This product is specially used for monitoring the biodecontamination process of vaporized hydrogen peroxide (VHP) and is suitable for the following scenarios:**

- △ Validation of biological decontamination effect in confined spaces such as isolators and sterile material transfer equipment
- △ Evaluation of disinfection effect of clean factory buildings, laboratories and other environments